



iSy® ODSecure® System –
reliable denture retention
made easy

For the reliable anchoring of dentures on implants


The iSy ODSecure System

Using the iSy ODSecure System (Over-Denture Secure System), dentures and partial dentures can be easily attached to dental implants. The **ODSecure Abutments** feature one of the flattest profiles in the industry, thus allowing for more flexible design of the denture base. They are

available in gingival heights ranging from 1.0 mm to 6.0 mm, which makes them a reliable fixation solution, even in challenging cases. And owing to the standardized iSy Platform, ODSecure Abutments also fit all iSy implant diameters. True to iSy.

Slim profile

Overall height of only 2.0 mm including housing cap and retention cap insert




Screwdrivers, hex

For tightening with a screwdriver, hex




Abrasion-resistant coating

Special titanium nitride coating for improved durability and strength



The **housing cap** corrects a divergence of up to 15°.



The **Xtend housing cap** corrects a divergence of up to 25°.



Product details

ODSecure Abutments

New

	Article	Art. No.	Ø	GH
	iSy® ODSecure® Abutment Material Titanium alloy	P2273.4010	4.0 mm	1.0 mm
		P2273.4020		2.0 mm
		P2273.4030		3.0 mm
		P2273.4040		4.0 mm
		P2273.4050		5.0 mm
		P2273.4060		6.0 mm

Accessories

	Article	Quantity	Color	Art. No.
	ODSecure® Attachment for Multi-unit abutment Ø 4.8 mm Material Titanium alloy	1	-	J2273.4801
	ODSecure® Assortment kit Material Titanium alloy/nylon/silicone	2	-	ODS-CAK
	ODSecure® Block-out spacer Material Silicone	4	clear	ODS-BS
	ODSecure® Analog Material Titanium alloy	2	-	ODSA
	ODSecure® Cap extractor tool Material stainless steel/nylon	1	-	ODSCT

The **ODSecure Assortment kit** includes numerous different components such as a housing cap with four retention cap inserts, one lab processing insert and a block-out spacer.



Accessories

	Article	Quantity	Color	Retention	Divergence	Art. No.
	ODSecure® Xtend housing cap incl. lab processing insert Material Titanium alloy/nylon	2	yellow anodized	-	0°-50°	ODS-XH
	ODSecure® housing cap incl. lab processing insert Material Titanium alloy/nylon	2	pink anodized	-	0°-30°	ODS-HCPM
	ODSecure® Retention cap insert Material Nylon	4	violet	HARD approx. 6 lbs / 2.7 kg	-	ODSRC-V
		4	clear	MEDIUM approx. 4 lbs / 1.8 kg	-	ODSRC-C
		4	pink	SOFT approx. 2.5 lbs / 1.2 kg	-	ODSRC-P
		4	yellow	EXTRA SOFT approx. 1.5 lbs / 0.6 kg	-	ODSRC-Y
	ODSecure® Lab processing insert Material Nylon	4	black	-	-	ODS-PM
	ODSecure® Impression coping Material Stainless steel/nylon	2	-	-	-	ODSIC
	ODSecure® Castable male burn-out Material Polystyrene	2	-	-	-	ODSCM

Small but extremely effective!

iSy Healing Caps and Impression Posts XS New

In future, the standardized iSy Platform offers even greater flexibility: In addition to the S, M, and L prosthetic profiles, the particularly slim XS profile can now also be used. This also makes the restoration of parti-

cularly narrow gaps with individual abutments possible, such as the use of the iSy ODSecure System for the restoration of edentulous and partially edentulous jaws with removable restorations.

	Article	Size	Art. No.	GH
	iSy® Healing caps, cylindrical sterile Material Titanium alloy	XS	P2015.0030	3.0 mm
			P2015.0045	4.5 mm
			P2015.0060	6.0 mm
	iSy® Impression posts, open tray incl. fixing screw Material Titanium alloy	XS	P2121.0001	-
	iSy® Impression posts, closed tray incl. fixing screw, impression cap and bite registration cap Material Titanium alloy/POM	XS	P2110.0001	-

BioHorizons Camlog Benelux B.V.

Pelmolenlaan 1 | 3447 GW Woerden | Nederland
info@camlog.nl | www.biohorizonscamlog.nl

Customer Service

Tel +31 (0)348 820 010 | info@camlog.nl

Headquarters

CAMLOG Biotechnologies GmbH | Margarethenstr. 38 | 4053 Basel | Switzerland
Phone +41 61 565 41 00 | Fax +41 61 565 41 01 | info@camlog.com | www.bhclgroup.com

Unless otherwise noted, all product names, whether printed in capital letters or bearing the trademark symbol, are trademarks of BioHorizons Inc. or CAMLOG Biotechnologies GmbH, its subsidiaries, affiliates, or its licensors. Manufacturer iSy® products: ALTATEC GmbH | Maybachstr. 5 | 71299 Wimsheim | Germany

Manufacturer of ODSecure® products:*

BioHorizons | 2300 Riverchase Center | Birmingham, AL 35244 | USA
ALTATEC GmbH | Maybachstr. 5 | 71299 Wimsheim | Germany

* This document refers to products from various approved manufacturers.
Please refer to the respective product label for the manufacturer's identification.

Valid from September 2025